# LuxaPrint Ortho Flex

DMG Digital Enterprises SE  
Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com  
092140/2021-01

## Gebrauchsinformation

Deutsch

### Produktbeschreibung

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EBPADMA, Urethanmethacrylat, TPO, Acrylat-Styren-Copolymer, Additive. Enthält TPO.

## Instructions for use

English

### Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP / SLA printers that work at wavelengths of 385 nm or 405 nm. Instruments with material parameters already defined in the system are recommended: e.g. DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ (Rapid Shape), P20+/P10+ capsule (Straumann). Select material variants according to the wavelength used by the printer (refer to »Packaging«).

### Intended purpose

Manufacturing of bleaching trays and transfer trays (in case of indirect bonding of brackets) using 3D printing.

### Limitations of use

Only insert LuxaPrint Ortho Flex intraorally in a fully polymerized state.

### Indications

- Bleaching trays: Tooth discoloration
- Transfer splints: Indirect bonding of brackets in orthodontic treatment.

### Contraindications

Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Intended users

Dentist, dental technician

### Notes for use

- During handling (e.g. transferring), the material should be exposed to daylight / ambient light only for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

### For optimum properties of the end product:

- Ensure that the instruments used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the instruments can cause defects on the printed object.
- Use the material parameters already defined in the system for the cleaning units DMG 3Dewash (DMG), RS wash (Rapid Shape) and P wash (Straumann).
- Allow objects to cure completely.
- Use the material parameters already defined in the system for the light-curing units DMG 3Decure (DMG), RS cure (Rapid Shape) and P cure (Straumann).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.
- Maximum curing depth during post-curing: For solid objects and when light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).

### Recommended use

**Note:** Observe the device and material parameters.

**Note:** Shake the bottle for at least 01:00 min before use.

**Note:** Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

**Note:** For optimal process reliability, scan the RFID code of the material (supported by DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ (Rapid Shape); P20+/P10+ capsule (Straumann).

## Finishing the printed object

<b>Production steps</b>	
Recommended drip-off time after raising the platform: 10:00 minutes.	
<b>Note:</b> Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform. When using the printers DMG 3Delite (DMG), D10+ capsule (Rapid Shape) and P10+ capsule (Straumann), start directly from step 2.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning instrument and specified program (e.g. DMG 3Dewash, RS wash or P wash).
2. (Alternatively)	Pre-clean the printed object with ethanol (approx. 96 %) or isopropanol (approx. 99 %) for max. 03:00 minutes in an ultrasonic cleaner (if necessary with a brush). Optional: Cleanse openings, holes and gap areas with compressed air. Rinse the printed object in a separate container with <b>fresh</b> ethanol (approx. 96 %) or isopropanol (approx. 99 %) for max. 02:00 minutes.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program (e.g. DMG 3Decure, RS cure or P cure).
4.	After removing the support, finish the surface.

## Interactions

Residues of methacrylate-based materials (uncured material residues) may influence the setting behavior of silicone impression and bite registration materials.

## Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

## Warnings / precautions

- Keep out of the reach of children!
- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid breathing in the fumes.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Observe the safety data sheet.
- Do not use the material where there are known allergies to any of the components or contact allergies.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Performance characteristics of the device

<b>Cured material</b>	
Shore-D hardness	≤ 80

<b>Uncured material</b>	
Viscosity	≤ 4.0 [Pa*s]
Light sensitivity to fluorescent lamps	≤ 05:00 [min]

## Storage / disposal

- Store in a dry place at room temperature (15 – 25 °C/59 – 77 °F) and protected from light
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

## Composition

EBPADMA, urethane methacrylate, TPO, styrene acrylate copolymer, additives. Contains TPO.

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